 **Lynx™** Telephone Transmitter


ENGLISH: Telephone transmitter for the Lynx system	page 3
DEUTSCH: Telefonsender für das Lynx System.....	Seite 8
FRANÇAIS: Émetteur téléphonique pour le système Lynx	page 13
ITALIANO: Trasmettitore telefonico per il sistema Lynx .	Pag. 18
NEDERLANDS: Telefoonzender voor het Lynx-systeem .	Pag. 23
NORSK: Telefonsender til Lynx-systemet	side 28
SVENSKA: Telefonsändare till Lynx-systemet.....	sid 33

Telephone transmitter for the Lynx system

The telephone transmitter forms part of Lynx – our wireless alert system for people with impaired hearing. Lynx uses various types of transmitters and receivers that alert you to noises such as a ringing telephone or doorbell, or the first signs of a fire.


The telephone transmitter sends signals wirelessly to the receiver in the system. One or more repeaters may also be installed in order to amplify the alarm signal and thus extend its range.


Connection examples


Select the appropriate connection. Not all parts on the diagrams are included in the package.

Example 1:

Example 2:**Wall mounting**

If you wish, you can mount the telephone transmitter on the wall.

**When the telephone rings**

When the telephone rings, the receiver installed in the Lynx system will be activated and will issue a warning by vibration and symbols, loud noise, or flashing light, depending on the type of receiver.

Read the user manual for the receiver installed!


Battery warning

The telephone transmitter is powered by a battery with a normal service life of approximately 3-5 years. When battery power begins to fail, the telephone transmitter will send out an alert. It sends a flat-battery warning along with the alarm to the receiver in the alert system.

Read the user manual for the receiver installed!

Changing the battery

1. Unscrew the two screws on the rear of the telephone transmitter, and release the rear cover by lifting it straight up.
2. Carefully push the old battery out of the holder.
3. Replace the battery with a new CR2032 Duracell type battery.
Note: The battery's positive terminal must be facing upwards.

4. Screw on the rear cover.

Changing the channel setting

The telephone transmitter sends signals wirelessly to the receiver installed in the Lynx system. It can transmit on a total of 15 different channels: 1-9, A-F. Channel 0 is reserved for a common fire alarm and should not be used.

Transmitters and receivers in the system must be set to the same channel so that they can communicate. They are all pre-set to channel 9, and then operate together without alteration.

In some cases, it may be necessary to change the channel setting. If two or more users live near one another, the alarms can interfere. This can be rectified by changing the channel settings so that they are different for everyone. Any channel can be selected except Channel 0.

If you need to change channels on the receiver, see the receiver user manual. If you need to change the channel on the telephone transmitter, proceed as follows:

The channel selector is on the rear. Use a small screwdriver and turn the yellow arrow towards the required channel.

Technical data

Modulation:	AM.
Frequency:	433,92 MHz.
Power supply:	3 V lithium battery (CR 2032). WARNING: Risk of explosion if fitted incorrectly!
Service life:	5 years (with 3 alarms per day). Maximum 70,000 alarms.
Range:	Up to 30 m indoors.
Channels:	16.
Connection:	Modular switch for adapter plug.
Dimensions:	60 x 40 x 23 mm.
Weight:	35 g.

CE Any issues relating to the EU Medical Device Directive 93/42/EEC should be directed to GN ReSound A/S.

Troubleshooting

Type of fault	Action
Telephone transmitter does not detect sound.	<ol style="list-style-type: none"> 1. Check that the telephone transmitter is connected correctly (see page 3). 2. Check that the battery is OK (see page 5).
The telephone transmitter gives a false alarm.	Someone in the vicinity has the same system installed. Change to another channel on all units.
Receiver fails to respond.	<ol style="list-style-type: none"> 1. Check the range. Try to reduce the distance between the telephone transmitter and the receiver or install a Repeater to extend the range. 2. Check that the receiver's battery is OK.

Cleaning and maintenance

Use a vacuum cleaner with a soft brush to clean the product on the outside. Use a damp cloth to wipe it on the outside. Never use cleaning products or solvents.

From time to time check that the cables and connections are intact.

The battery/batteries has/have a normal operating life of about 3-5 year(s). Refer to page 5 for how to change the battery.

Repairs

If your product malfunctions, it must be repaired by a qualified technician. Do not attempt to open the case of the device since this would invalidate the warranty. If your product requires service, please contact your hearing care professional for assistance.

Waste disposal

This product contains electrical or electronic equipment and should be disposed of carefully in the interest of your safety and the environment. Please contact your local hearing care professional concerning disposal of the product.

General precautions

- Low-energy radiation (eg from burglary alarms, room surveillance systems, radio equipment, mobile telephones) will not damage the device. They could, however, momentarily affect the sound quality or create noise in the device.
- Warning: External devices connected to the product must be safe according to the requirements of EN 60601-1-1, EN 60335-1 or EN 60065, as appropriate.

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GN Resound AB